


Product Name: Recombinant Human CALD1 Protein
Catalog #: PEH7658



Summary

Name	CALD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05682
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KREERRKVL EEEQRRKQEE ADRKLREEEE KRRLKEEIER RRAEAAEKRQ KMPEDGLSDD KKPFCFTPK GSSLKIEERA EFLNKSQKS SGVKSTHQAA IVSKIDSRLE QYTS AIEGTK SAKPTKPAAS DLPVPAEGVR NIKSMWEKGN VFSSPTAAGT PNKETAGLKV GVSSRINWL TKTPDGNKSP APKPSDLRPG DVSSKRNLWE KQSVDKVTSP

image



Background

Alternative Names Caldesmon (CDM)

Note

Product Name: Recombinant Human CALD1 Protein
Catalog #: PEH7658



For Research Use Only , Not for Diagnostic Use