


Product Name: Recombinant Human RNASEL Protein
Catalog #: PEH7660



Summary

Name	RN5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05823
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEVLKILLDE MGADVACDN MGRNALIHAL LSSDDSDVEA ITHLLLDHGA DVNVRGERGK TPLILAVEKK HLGLVQRLLE QEHIENDTD

image



Background

Alternative Names	2-5A-dependent ribonuclease (2-5A-dependent RNase) (EC 3.1.26.-) (Ribonuclease 4) (Ribonuclease L) (RNase L)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use