


Product Name: Recombinant Human TJP1 Protein
Catalog #: PEH7661



Summary

Name	ZO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q07157
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAMEETAIWE QHTVTLHRAP GFGFGIAISG GRDNPHFQSG ETSIVISDVL KGGPAEGQLQ ENDRVAMVNG VSMDNVEHAF AVQQLRKSGK NAKITIRKK KVQIPVSRPD PEPVSDNEED SYDEEIHDPGR SGRSGVVNRR SEKIWPRDRS ASRERSLSPR SDRRSVASSQ PAKPTKVTLV KSRKNEEYGL RLASHIFVKE ISQDSLAARD GNIQEGDVVL KINGVTENM SLTDAKTIE RSKGKLKMVV

image



Background

Alternative Names	Tight junction protein ZO-1 (Tight junction protein 1) (Zona occludens protein 1) (Zonula occludens protein 1)
--------------------------	--

Note

Product Name: Recombinant Human TJP1 Protein
Catalog #: PEH7661



For Research Use Only , Not for Diagnostic Use