


Product Name: Recombinant Human ITK Protein
Catalog #: PEH7662



Summary

Name	ITK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08881
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NASKKPLPPT PEDNRRPLWE PEETVVIALY DYQTNDPQEL ALRRNEEYCL LDSSEIHWWR VQDRNGHEGY VPSSYLVEKS PNNLETYEWY NKSISRDKAE KLLLDTGKEG AFMVRDSRTA GTYTVSVFTK AVVSENNPCI KHYHIKETND NPKRYYYAEK YVFDSIPLLI NYHQHNGGGL

image



Background

Alternative Names	Tyrosine-protein kinase ITK/TSK (EC 2.7.10.2) (Interleukin-2-inducible T-cell kinase) (IL-2-inducible T-cell kinase) (Kinase EMT) (T-cell-specific kinase) (Tyrosine-protein kinase Lyk)
--------------------------	--

Note

Product Name: Recombinant Human ITK Protein
Catalog #: PEH7662



For Research Use Only , Not for Diagnostic Use