


Product Name: Recombinant Human DMPK Protein
Catalog #: PEH7663



Summary

Name	DMPK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q09013
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MELEAEQ LLE PHVQAPSLEP SVSPQDETAE VAVPAAVPAA EAEAEVTLRE LQEAL EEEVL TRQSL SREME AIRTDNQNFA SQLREAEARN RDLEAHVRQL QERMELLQAE

image



Background

Alternative Names Myotonin-protein kinase (MT-PK) (EC 2.7.11.1) (DM-kinase) (DMK) (DM1 protein kinase) (DMPK) (Myotonic dystrophy protein kinase)

Note

For Research Use Only , Not for Diagnostic Use