


Product Name: Recombinant Human NEXN Protein
Catalog #: PEH7664



Summary

Name	NEXN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0ZGT2
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLSSSKPVPK TYVPKLGKGD VKDKFEAMQR AREERNQRRS RDEKQRRKEQ YIREREWNRR KQEIKEMLAS DDEEDVSSKV EKAYVPKLTG TVKGRFAEME KQRQEEQRKR TEEERKRRIE QDMLEKRRIQ RELAKRAEQI EDINNTGTES ASEEGDSSL ITVVPVKS YK TSGKMKNFE DLEKEREKE RIKYEEDKRI RYEEQRPSLK EAKCLSLVMD DEIESEAKKE SLSPGKLKLT FEELERQRQE NRKKQAEAEA RKRLEEEKRA FEEARRQMVN EDEENQDTAK IFKGYRPGKL KLSFEEMERQ RREDEKRKAE EEARRRIEEE KKAFAEARRN MVVDDDSPEN YKTISQEFLT

image



Background

Alternative Names	Nexilin (F-actin-binding protein) (Nelin)
--------------------------	---

Product Name: Recombinant Human NEXN Protein
Catalog #: PEH7664



Note

For Research Use Only , Not for Diagnostic Use