



Product Name: Recombinant Human NFX1 Protein
Catalog #: PEH7666



Summary

Name	NFX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12986
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FKFNTDAAEF IPQEKKNSSL NCGTQRRLLS NRIGRRNYSS PPPCHLSRQV PYDEISAVHQ HSYHPSGSKP KSQQTSTFQSS PCNKSPKSHG LQNQPWQKLR NEKHHRVKK AQSLAEQTSN TAGLESSTRS ESGTDLREHS PSESEKEVVG ADPRGAKPKK ATQFVYSYGR GPKVKGKLLC EWSNRTPPKP EDAGPESTKP VGVFHPDSSE ASSRKGVLDT YGARRNEQRR YPQKRPPWEV EGARPRPGRN PPKQEGHRHT NAGHRNNMGP IPKDDLNERP AKSTCDSENH AVINKSSRRV DQEKCTVRRQ DPQVVSPFSR GKQNHVLKNV ETHTGSLIEQ

image



Background

Alternative Names	Transcriptional repressor NF-X1 (EC 2.3.2.-) (Nuclear transcription factor, X box-binding protein 1)
--------------------------	--

Product Name: Recombinant Human NFX1 Protein
Catalog #: PEH7666



Note

For Research Use Only , Not for Diagnostic Use