



Product Name: Recombinant Human FLII Protein
Catalog #: PEH7667



Summary

Name	FLII_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13045
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WVPVEYEEEE KKEDKEEKAE GKEGEEATAE AEEKQPEEDF QCIVYFWQGR EASNMGWLTF TFSLQKKFES LFP GKLEVVR MTQQQENPKF LSHFKRKFII HRGKRKAVQG AQQPSLYQIR TNGSALCTRC IQINTDSSLL NSEFCFILKV PFESDNQGI VYAWVGRASD PDEAKLAEDI LNTMFDTSYS KQVINEGEEP ENFFWVGIGA QKPYDDDAEY MKHTRLFRCS NEKGYFAVTE KCSDFCQDDL ADDDIMLLDN GQEVYMWVGT QTSQVEIKLS LKACQVYIQH MRSKEHERPR RLRLVRKGN E QHAFTRCFHA

image



Background

Alternative Names Protein flightless-1 homolog

Product Name: Recombinant Human FLII Protein
Catalog #: PEH7667



Note

For Research Use Only , Not for Diagnostic Use