


Product Name: Recombinant Human COASY Protein
Catalog #: PEH7668



Summary

Name	COASY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13057
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLVVGVADKD LLKSKLLPEL LQPYTERVEH LSEFLVDIKP SLTFDVIPLL DPYGPAGSDP SLEFLVSEE TYRGGMAINR FRLNDLEEL ALYQIQLLKD LRHTENEEDK

image



Background

Alternative Names	Bifunctional coenzyme A synthase (CoA synthase) (NBP) (POV-2) [Includes: Phosphopantetheine adenyltransferase (EC 2.7.7.3) (Dephospho-CoA pyrophosphorylase) (Pantetheine-phosphate adenyltransferase) (PPAT); Dephospho-CoA kinase (DPCK) (EC 2.7.1.24) (Dephosphocoenzyme A kinase) (DPCOAK)]
--------------------------	---

Product Name: Recombinant Human COASY Protein
Catalog #: PEH7668



Note

For Research Use Only , Not for Diagnostic Use