



Product Name: Recombinant Human TBX2 Protein
Catalog #: PEH7672



Summary

Name	TBX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13207
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKELWDQFHK LGTEMVITKS GRRMFPPFKV RVSGLDKKAK YILLMDIVAA DDCRYKFHNS RWMVAGKADP EMPKRMVYIHP DSPATGEQWM AKPVAFHKLK LTNNISDKHG FTILNSMHKY QPRFHIVRAN DILKLPYSTF RTYVFPETDF IAVTAYQNDK ITQLKIDNNP FAKGFRDTGN GRREKRKQLT LPSLRLYEEH CKPERDGAES DASSCDPPPA REPPTSPGAA PSPLRLHRAR AEEKSCAADS DPEPERLSEE RAGAPLGRSP

image



Background

Alternative Names T-box transcription factor TBX2 (T-box protein 2)

Product Name: Recombinant Human TBX2 Protein
Catalog #: PEH7672



Note

For Research Use Only , Not for Diagnostic Use