


Product Name: Recombinant Human MAP3K1 Protein
Catalog #: PEH7674



Summary

Name	M3K1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13233
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLPSSNIHRP KPSRTPGNT SKQGDPKNS MTLDLNSSSK CDDSGCASN SSNAVIPSDE TVFTPVEEK RLDVNTELNS SIEDLLEASM PSSDTTVTFK SEVAVLSPEK AENDDTYKDD VNHNQKCKEK MEAEELALAL IAMAMSASQD ALPIVPQLQV ENGEDIILIQ QDTPETLPGH

image



Background

Alternative Names	Mitogen-activated protein kinase kinase kinase 1 (EC 2.7.11.25) (MAPK/ERK kinase kinase 1) (MEK kinase 1) (MEKK 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MAP3K1 Protein
Catalog #: PEH7674

