

Product Name: Recombinant Human PABPC4 Protein
Catalog #: PEH7676



Summary

Name	PABP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13310
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	<pre>EVDDESLKEL FSQFGKTLV KVMRDPNGKS KGFGFVSYEK HEDANKAVEE MNGKEISGKI IFVGRAQKKV ERQAEKRRKF EQLKQERISR YQGVNLYIKN LDDTIDDEKL RKEFSPFGSI TSAKVMLEDG RSKGFGFVCF SSPEEATKAV TEMNGRIVGS KPLYVALAQR KEERKAHLTN QYMQRVAGMR ALPANAILNQ FQPAAGGYFV PAVPQAQGRP PYYTPNQLAQ MRPNPRWQQG GRPQGFQGMF SAIRQSGPRP TLRHLAPTGS ECPDRLAMDF GGAGAAQQGL TDSCQSGGVP</pre>

image

SDS

Background

Alternative Names Polyadenylate-binding protein 4 (PABP-4) (Poly(A)-binding protein 4) (Activated-platelet protein 1) (APP-1) (Inducible poly(A)-binding protein) (iPABP)

Note

For Research Use Only , Not for Diagnostic Use