


Product Name: Recombinant Human EIF3I Protein
Catalog #: PEH7678



Summary

Name	EIF3I_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13347
Host	E.coli
Species	Human
Molecular Mass	79 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKPILLQGHE RSITQIKYNR EGDLLFTVAK DPIVNVWYSV NGERLGTYMG HTGAVWCVDA DWDTKHVLTG SADNSCRLWD CETGKQLALL KTNSAVRTCG FDFGGNIIMF STDKQMGYQC FVSFFDLRDP SQIDNNEPYM KIPCNDKIT SAVWGPLGEC IIAGHESGEL NQYSAKSGEV LVNVKEHSRQ INDIQLSRDM TMFVTASKDN TAKLFDSTTL EHQTFRTER PVNSAALSPN YDHVVLGGGQ EAMDVTTTST RIGKFEARFF HLAFAEEFGR VKGHFGPINS VAFHPDGKSY SSGGEDGYVR IHYFDPQYFE FEFEA

image



Background

Alternative Names	Eukaryotic translation initiation factor 3 subunit I (eIF3i) (Eukaryotic translation initiation factor 3 subunit 2) (TGF-beta receptor-interacting
--------------------------	--

Product Name: Recombinant Human EIF3I Protein
Catalog #: PEH7678



protein 1) (TRIP-1) (eIF-3-beta) (eIF3 p36)

Note

For Research Use Only , Not for Diagnostic Use