


Product Name: Recombinant Human EIF3I Protein
Catalog #: PEH7679



Summary

Name	EIF3I_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13347
Host	E.coli
Species	Human
Molecular Mass	70 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DWDTKHVLTG SADNSCRLWD CETGKQLALL KTNSAVRTCG FDFGGNIIMF STDKQMGYQC FVSFFDLRDP SQIDNNEPYM KIPCNSKIT SAVWGPLGEC IIAGHESGEL NQYSAKSGEV LVNVKEHSRQ INDIQLSRDM TMFVTASKDN TAKLFDSTTL EHQKTFRTER PVNSAALSPN YDHVVLGGGQ EAMDVTTTST RIGKFEARFF HLAFAEEEFGR VKGHFGPINS VAFHPDGKSY

image



Background

Alternative Names	Eukaryotic translation initiation factor 3 subunit I (eIF3i) (Eukaryotic translation initiation factor 3 subunit 2) (TGF-beta receptor-interacting protein 1) (TRIP-1) (eIF-3-beta) (eIF3 p36)
--------------------------	--

Product Name: Recombinant Human EIF3I Protein
Catalog #: PEH7679



Note

For Research Use Only , Not for Diagnostic Use