


Product Name: Recombinant Human TNK1 Protein
Catalog #: PEH7681



Summary

Name	TNK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13470
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPLERMKGIS RSLESVLSLG PRPTGGGSSP PEIRQARAVP QGPPGLPPRP PLSSSSPQPS QPSRERLPWP KRKPPHNHPM

image



Background

Alternative Names	Non-receptor tyrosine-protein kinase TNK1 (EC 2.7.10.2) (CD38 negative kinase 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use