



Product Name: Recombinant Human RIPK1 Protein
Catalog #: PEH7684



Summary

Name	RIPK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q13546
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KMKSSDFLES AELDSGGFGK VSLCFHRTQG LMIMKTVYKG PNCIEHNEAL LEEAKMMNRL RHSRVVKLLG VIIEEGKYSL VMEYMEKGNL MHVLKAEMST PLSVKGRIIL EIIEGMCYLH GKGVIHKDLK PENILVDNDF HIKIADLGLA SFKMWSKLNN EEHNELREVD GTAKKNGGTL YYMAPEHLND VNAKPTEKSD VYSFAVVLWA IFANKEPYEN AICEQQLIMC IKSGNRPDVD DITEYCPREI ISLMKLCWEA NPEARPTFPG IEEKFRPFYL

image



Background

Alternative Names	Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP) (Receptor-interacting protein 1) (RIP-1)
--------------------------	--

Product Name: Recombinant Human RIPK1 Protein
Catalog #: PEH7684



Note

For Research Use Only , Not for Diagnostic Use