


Product Name: Recombinant Human SNTB1 Protein
Catalog #: PEH7692



Summary

Name	SNTB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13884
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRSGLLLEVLV RDRWHKVLVN LSEDALVLSS EEGAAAYNGI GTATNGSFRCR GAGAGHPGAG GAQPPDSPAG VRTAFTDLPE QVPESISNQK RGVKVLKQEL GGLGISIKGG KENKMPILIS KIFKGLAADQ TQALYVGDAI LSVNGADLRD ATHDEAVQAL KRAGKEVLLE VKYMREATPY VKKGSPVSEI GWETPPPEPSP RLGGSTSDPP SSQSFSFHRD RKSIPKMCY VTRSMALADP ENRQLEIHSP DAKHTVILRS KDSATAQAWF

image



Background

Alternative Names	Beta-1-syntrophin (59 kDa dystrophin-associated protein A1 basic component 1) (DAPA1B) (BSYN2) (Syntrophin-2) (Tax interaction protein 43) (TIP-43)
--------------------------	---

Product Name: Recombinant Human SNTB1 Protein
Catalog #: PEH7692



Note

For Research Use Only , Not for Diagnostic Use