


Product Name: Recombinant Human CNGB1 Protein
Catalog #: PEH7694



Summary

Name	CNGB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14028
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPPSPAKSDT LIVPSSASGT HRKKLPSEDD EAEELKALSP AESPVVAWSD PTTPKDTDGQ DRAASTASTN SAIINDRLQE LVKLFKERTE KVKEKLIDPD VTSDEESPKP SPAKKAPEPA PDTKPAAEAP VEEHYCDML CCKFKHRPWK

image



Background

Alternative Names	Cyclic nucleotide-gated cation channel beta-1 (Cyclic nucleotide-gated cation channel 4) (CNG channel 4) (CNG-4) (CNG4) (Cyclic nucleotide-gated cation channel gamma) (Cyclic nucleotide-gated cation channel modulatory subunit) (Cyclic nucleotide-gated channel beta-1) (CNG channel beta-1) (Glutamic acid-rich protein) (GARP)
--------------------------	--

Product Name: Recombinant Human CNGB1 Protein
Catalog #: PEH7694



Note

For Research Use Only , Not for Diagnostic Use