


Product Name: Recombinant Human TTL12 Protein
Catalog #: PEH7695



Summary

Name	TTL12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14166
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAERSSPGQT PEEGAQALAE FAALHGPAIR ASGVPERYWG RLLHKLEHEV FDAGEVFGIM QVEEVEEED EAAREVRKQQ PNPGNELCYK VIVTRESGLQ AAHPNSIFLI DHAWTCRVEH ARQQLQQVPG LLHRMANLMG IEFHGELPST EVALVLEEM WKFNQTYQLA HGTAEEKMPV WYIMDEFGSR IQHADVPSFA TAPFFYMPQQ VAYTLLWPLR DLDTGEEVTR DFAYGETDPL IRKCMLLPWA PTDMLDLSSC TPEPPAEHYQ AILEENKEKL

image



Background

Alternative Names	Tubulin--tyrosine ligase-like protein 12 (Inactive tubulin--tyrosine ligase-like protein 12)
--------------------------	--

Product Name: Recombinant Human TTLL12 Protein
Catalog #: PEH7695



Note

For Research Use Only , Not for Diagnostic Use