

Product Name: Recombinant Human TFDP1 Protein
Catalog #: PEH7696



Summary

Name	TFDP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14186
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FEYLFNFDNT FEIHDDIEVL KRMGMACGLE SGSCSAEDLK MARSLVPKAL EPYVTEMAQG TVGGVFITTA GSTSNGTRFS ASDLTNGADG MLATSSNGSQ YSGSRVETPV SYVGEDDEED

image

SDS

Background

Alternative Names	Transcription factor Dp-1 (DRTF1-polypeptide 1) (DRTF1) (E2F dimerization partner 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use