


Product Name: Recombinant Human KPNB1 Protein
Catalog #: PEH7701



Summary

Name	IMB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14974
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSPDRLELEA AQKFLERA AV ENLPTFLVEL SRVLANPGNS QVARVAAGLQ IKNSLTSKDP DIKAQYQQRW LAIDANARRE VKNYVLQTLG TETYPSSAS QCVAGIACAE IPVNQWP ELI PQLVANVTNP NSTEHMKEST LEAIGYICQD IDPEQLQDKS NEILTAIIQG MRKEEPSNNV KLAATNALLN SLEFTKANFD KESERHFIMQ VVCEATQCPD TRVRVAALQN LVKIMSLYYQ YMETYMG PAL FAITIEAMKS DIDEVALQGI EFWSNVCDEE MDLAI EASEA AEQGRPPEHT SKFYAKGALQ YLVPILTQTL TKQDENDDDD DWNPCKAAGV CLMLLATCCE DDIVPHVLPF IKEHIKNPDW RYRDAAVMAF GCILEGPEPS QLKPLVIQAM PTLIELMKDP SVVVRDTAAW TVGRICELLP

image



Background

Product Name: Recombinant Human KPNB1 Protein
Catalog #: PEH7701



Alternative Names

Importin subunit beta-1 (Importin-90) (Karyopherin subunit beta-1) (Nuclear factor p97) (Pore targeting complex 97 kDa subunit) (PTAC97)

Note

For Research Use Only , Not for Diagnostic Use