


Product Name: Recombinant Human NR1I3 Protein
Catalog #: PEH7702



Summary

Name	NR1I3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14994
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CVVCGDQATG YHFNALTCEG CKGFFRRTVS KSIPTCPFA GSCEVSKTQR RHCPACRLQK CLDAGMRKDM ILSAEALALR RAKQAQRRQA

image



Background

Alternative Names	Nuclear receptor subfamily 1 group I member 3 (Constitutive activator of retinoid response) (Constitutive active response) (Constitutive androstane receptor) (CAR) (Orphan nuclear receptor MB67)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use