


Product Name: Recombinant Human GAPVD1 Protein
Catalog #: PEH7703



Summary

Name	GAPD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14C86
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSQLITLVNF MKSVMSGDQL REDRMALDNL LANLPPAKPG KSSSLEMTPY NTPQLSPATT PANKKNRLPI ATRSRRTNM LMDLHMDHEG SSQETIQEVQ PEEVLVISLG TGPQLTPGMM SENEVLNMQ L SDGGQGDVPV DENKLHGKPD KTLRFSLCSD NLEGISEGPS NRSNSVSSLD LEGESVSELG AGPSGSNGVE ALQLEHEQA TTQDNLDDKL RKFEIRDMMG LTDDRDISET VSETWSTDVL GSDFDPNIDE DRLQEIAGAA

image



Background

Alternative Names	GTPase-activating protein and VPS9 domain-containing protein 1 (GAPex-5) (Rab5-activating protein 6)
--------------------------	---

Product Name: Recombinant Human GAPVD1 Protein
Catalog #: PEH7703



Note

For Research Use Only , Not for Diagnostic Use