

Product Name: Recombinant Human TNIP1 Protein
Catalog #: PEH7704



Summary

Name	TNIP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15025
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPGGSVPSGE ASAAFERLVK ENSRLKEKMQ GIKMLGELLE ESQMEATRLR QKAEELVKDN ELLPPSPSL GSFDPLAELT GKDSNVTASP TAPACPSDKP APVQKPPSSG TSSEFEVVTP EEQNSPESSS HANAMALGPL PREDGNLMLH LQRLETTLSV CAEEDPDHGQL FTHLGRMALE FNRLASKVHK NEQRTSILQT LCEQLRKENE ALKAKLDKGL EQRDQAAERL REENLELKKL LMSNGNKEGA SGRPGSPKME GTGKKAVAGQ

image



Background

Alternative Names	TNFAIP3-interacting protein 1 (A20-binding inhibitor of NF-kappa-B activation 1) (ABIN-1) (HIV-1 Nef-interacting protein) (Nef-associated factor 1) (Naf1) (Nip40-1) (Virion-associated nuclear shuttling protein) (VAN)
--------------------------	--

Product Name: Recombinant Human TNIP1 Protein
Catalog #: PEH7704



(hVAN)

Note

For Research Use Only , Not for Diagnostic Use