


Product Name: Recombinant Human EIF4H Protein
Catalog #: PEH7705



Summary

Name	IF4H_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15056
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MADFDTYDDR AYSSFGGGRG SRGSAGGHGS RSQKELPTEP PYTAYVGNLP FNTVQGDIDA IFKDLIRSIV RLVRDKDSDK FKGFCYVEFD EVDSLKEALT YDGALLGDRS LRVDIAEGRK QDKGGFGFRK GGPDDRGMGS SRESRGGWDS RDDFNSGFRD DFLGGRGGSR PGDRRTGPPM GSRFRDGPPL RGSNMDFREP TEEERAQRPR LQLKPRTVAT PLNQVANPNS AIFGGARPRE EVVQKEQE

image



Background

Alternative Names	Eukaryotic translation initiation factor 4H (eIF-4H) (Williams-Beuren syndrome chromosomal region 1 protein)
--------------------------	--

Note

Product Name: Recombinant Human EIF4H Protein
Catalog #: PEH7705



For Research Use Only , Not for Diagnostic Use