



Product Name: Recombinant Human PDIA6 Protein
Catalog #: PEH7710



Summary

Name	PDIA6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15084
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LYSSDDVIE LTPSNFNREV IQSDSLWLVE FYAPWCGHCQ RLTPWKKAA TALKDVVKVG AVDADKHSL GGQYGVQGFP TIKIFGSNKN RPEDYQGGRT GEAIVDAALS ALRQLVKDRL GGRSGGYSSG KQGRSDSSK KDVIELTDDS FDKNVLDSFD VWMVEFYAPW CGHCKNLEPE WAAAASEVKE QTKGKVKLAA VDATVNQVLA SRYGIRGFPT IKIFQKGESP VDYDGGRTS DIVSRALDLF SDNAPPELL EIINEDIAKR TCEEHQLCVV AVLPHILDTG AAGRNSYLEV LLKLADKYKK KMWGWLWTEA GAQSELETAL GIGGFGYPAM AAINARKMKF ALLKGSFSEQ GINEFLRELS FGRGSTAPVG GGAFTPTVER EPWDGRDGEL PVEDDIDLSD VELDDLKDE L

image



Background

Product Name: Recombinant Human PDIA6 Protein
Catalog #: PEH7710



Alternative Names

Protein disulfide-isomerase A6 (EC 5.3.4.1) (Endoplasmic reticulum protein 5) (ER protein 5) (ERp5) (Protein disulfide isomerase P5) (Thioredoxin domain-containing protein 7)

Note

For Research Use Only , Not for Diagnostic Use