



Product Name: Recombinant Human PTPA Protein
Catalog #: PEH7711



Summary

Name	PTPA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15257
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAEGERQPPP DSSEEAPPAT QNFIIPKKEI HTVPDMGKWK RSQAYADYIG FILTLNEGVK GKKLTFEYRV SEMWNEVHEE KEQAAKQSVS CDECIPLPRA GHCAPSEAIE KLVALLNTLD RWIDETPPVD QPSRFGNKAY RTWYAKLDEE AENLVATVVP THLAAVPEV AVYLKESVGN STRIDYGTGH EAAFAAFLCC LCKIGVLRVD DQIAIVFKVF NRYLEV MRKL QKTYRMEPAG SQGVWGLDDF QFLPFIWGSS QLIDHPYLEP RHFVDEKAVN ENHKDYMFLC CILFITEMKT GPFAEHSNQL WNISAVPSWS KVNQGLIRMY KAECLEKFPV IQHFKFGSLL PIHPVTSG

image



Background

Alternative Names	Serine/threonine-protein phosphatase 2A activator (EC 5.2.1.8) (PP2A,
--------------------------	---

Product Name: Recombinant Human PTPA Protein
Catalog #: PEH7711



subunit B', PR53 isoform) (Phosphotyrosyl phosphatase activator) (PTPA)
(Serine/threonine-protein phosphatase 2A regulatory subunit 4)
(Serine/threonine-protein phosphatase 2A regulatory subunit B')

Note

For Research Use Only , Not for Diagnostic Use