


Summary

Name	PCBP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15366
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NVTLTIRLLM HGKEVGSIIIG KKGESVKKMR EESGARINIS EGNCPERIIT LAGPTNAIFK AFAMIIDKLE EDISSMTNS TAASRPPVTL RLVVPASQCG SLIGKGGCKI KEIRESTGAQ VQVAGDMLPN STERAITIAG IPQSIIECVK QICVVMLETL SQSPPKGVTI PYRPKPSSSP VIFAGGQDRY STGSDSASFP HTTSPMCLNP DLEGPPELAY TIQGQYAIQ PDLTKLHQLA MQQSHFPMTH GNTGFSGIES SSPEVKGYWG LDASAQTTSHT ELTIPNDLIG CIIGRQGAKI NEIRQMGAQ IKIANPVEGS TDRQVTITGS AASISLAQYL INVRLSSETG

image



Background

Alternative Names	Poly(rC)-binding protein 2 (Alpha-CP2) (Heterogeneous nuclear ribonucleoprotein E2) (hnRNP E2)
--------------------------	--

Product Name: Recombinant Human PCBP2 Protein
Catalog #: PEH7713



Note

For Research Use Only , Not for Diagnostic Use