


Product Name: Recombinant Human SKIC2 Protein
Catalog #: PEH7714



Summary

Name	SKIV2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15477
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEFPSRKDSK AHEQALAEIT KRLGALEEPD MTGQLVDLPE YYSWGEELTE TQHMIQRRIM ESVNGLKSL AGRVVVKNQ EHHNALGVIL QVSSNSTSRV FTTLVLCDKP LSQDPQDRGP ATA EVPYPDD LVGFKLFLPE GPCDHTVVKL QPGDMAAITT KVLRVNGEKI LEDFSKRQQP KFKKDPPLAA VTTAVQELLR LAQAHPAGPP TLDPVNDLQL KDMSVVEGGL RARKLEELIQ GAQC VHSRPF PAQYLKLRER MQIQKEMERL

image



Background

Alternative Names Helicase SKI2W (Ski2) (EC 3.6.4.-) (Helicase-like protein) (HLP)

Product Name: Recombinant Human SKIC2 Protein
Catalog #: PEH7714



Note

For Research Use Only , Not for Diagnostic Use