


Product Name: Recombinant Human PTPN14 Protein
Catalog #: PEH7718



Summary

Name	PTN14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15678
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	REENRVDATR VPMDERFTL KKKLEEGMVF TEYEQIPKKK ANGIFSTAAL PENAERSRIR EVVPYEENRV ELIPTKENNT GYINASHIKV VVGGAEWHYI ATQGPLPHTC HDFWQMVWEQ GVNVIAMVTA EEEGGRKSH RYWPKLGSKH SSATYGKFKV TTKFRD SVC YATTGLKVKH LLSGQERTVW HLQYTDWPDH GCPEDVQGFL SYLEEQSVR RHTNSMLEGT KNRHPPIVVH CSAGVGRTGV LILSELMIYC LEHNEKVEVP MMLRLLREQR MFMIQTIAQY KFYVQVLIQF

image



Background

Alternative Names	Tyrosine-protein phosphatase non-receptor type 14 (EC 3.1.3.48) (Protein-tyrosine phosphatase pez)
--------------------------	--

Product Name: Recombinant Human PTPN14 Protein
Catalog #: PEH7718



Note

For Research Use Only , Not for Diagnostic Use