


Product Name: Recombinant Human ITSN1 Protein
Catalog #: PEH7721



Summary

Name	ITSN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15811
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLEKKLPVTF EDKKRENFER GNLELEKRRQ ALLEQQRKEQ ERLAQLERAE QERKERERQE QERKRQLELE KQLEKQRELE RQREEERRKE IERREAAKRE LERQRQLEWE RNRRQELLNQ RNKEQEDIVV LKAKKKTLEF ELEALNDKKH QLEGKLDIR CRLTTQRQEI ESTNKSRELR IAEITHLQQQ

image



Background

Alternative Names Intersectin-1 (SH3 domain-containing protein 1A) (SH3P17)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ITSN1 Protein
Catalog #: PEH7721

