


Product Name: Recombinant Human USF2 Protein
Catalog #: PEH7723



Summary

Name	USF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15853
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRDERRRAQH NEVERRRRDK INNWIVQLSK IIPDCNADNS KTGASKGGIL SKACDYIREL RQTNQRMQET FKEAERLQMD NELLRQQIEE LKNENALLRA QLQQHNLEMV

image



Background

Alternative Names	Upstream stimulatory factor 2 (Class B basic helix-loop-helix protein 12) (bHLHb12) (FOS-interacting protein) (FIP) (Major late transcription factor 2) (Upstream transcription factor 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USF2 Protein
Catalog #: PEH7723

