


Product Name: Recombinant Human OCLN Protein
Catalog #: PEH7727



Summary

Name	OCLN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16625
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDRYDKSNIL WDKEHIYDEQ PPNVEEWVKN VSAGTQDVPS PPSDYVERVD SPMAYSSNGK VNDKRFYPES SYKSTPVPEV VQELPLTSPV DDFRQPRYSS GGNFETPSKR APAKGRAGRS KRTEQDHYET DYTGGESCD ELEEDWIREY PPITSDQQRQ LYKRNFD TGL QEYKSLQSEL DEINKELSRL DKELDDYREE SEEYMAAADE YNRLKQVKGS ADYKSKKNHC KQLKSKLSHI KKMVG DYDRQ

image



Background

Alternative Names Occludin

Note

Product Name: Recombinant Human OCLN Protein
Catalog #: PEH7727



For Research Use Only , Not for Diagnostic Use