


Product Name: Recombinant Human GUK1 Protein
Catalog #: PEH7728



Summary

Name	KGUA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16774
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSGPRPVVLS GPSGAGKSTL LKRLLEHSG IFGFSVSHTT RNPRPGEENG KDYYFVTREV MQRDIAAGDF IEHAEFSGNL YGTSKVAVQA VQAMNRCVL DVDLQGVRNI KATDLRPIYI SVQPPSLHVL EQRLRQRNTE TEESLVKRLA AAQADMESK EPGLFDVVII NDSLQAYAE LKEALSEEIK KAQRTGA

image



Background

Alternative Names Guanylate kinase (EC 2.7.4.8) (GMP kinase) (Guanylate kinase 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GUK1 Protein
Catalog #: PEH7728

