



Product Name: Recombinant Human UGP2 Protein
Catalog #: PEH7730



Summary

Name	UGPA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16851
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AMSQDGASQF QEVIRQELEL SVKKELEKIL TTASSHEFEH TKKDLGFRK LFHRFLQKEG PSVDWGGKIQR PPEDSIQPYE

image



Background

Alternative Names	UTP--glucose-1-phosphate uridylyltransferase (EC 2.7.7.9) (UDP-glucose pyrophosphorylase) (UDPGP) (UGPase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use