

Product Name: Recombinant Human INF2 Protein
Catalog #: PEH7732



Summary

Name	INF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q27J81
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEQLEEGGPR PLERRSSWYV DASDVLTTED PQCPQPLEGA WPVTLGDAQA LKPLKFSSNQ PPAAGSSRQD AKDPTSLLGV LQAEADSTSE GLEDAVHSRG ARPPAAGPGG DEDEDEEDTA PESALDTSLD KSFSEDAVTD SSGSGTLPRA RGRASKGTGK RRRKRPSRSQ EEVPPDSDDN

image

SDS

Background

Alternative Names Inverted formin-2 (HBEBP2-binding protein C)

Note

For Research Use Only , Not for Diagnostic Use