


**Product Name: Recombinant Human GON4L Protein**  
**Catalog #: PEH7735**



## Summary

<b>Name</b>	GON4L_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q3T8J9
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	30 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DGPPPHSPET P QFPPTTGAV LYTVKRNQVG PEVRSCPAS PRLQKEREQQ KAVSESEALM LVWDASETEK LPGTVEPPAS FLSPVSSKTR DAGRRHVSGK PDTQERWLPS SRARVKTRDR TCPVHESPSG IDTSETSPKA PRGGLAKDSG TQAKGPEGEQ QPKAAEATVC ANNSKVSSTG EKVVLWTREA DRVILTMCQE QGAQPQTFNI ISQQLGNKTP AEVSHRFREL MQLFHTACEA SSEDEDDATS TSNADQLSDH GDLLSEEELD

## image



## Background

**Alternative Names** GON-4-like protein (GON-4 homolog)

**Product Name: Recombinant Human GON4L Protein**  
**Catalog #: PEH7735**

---



**Note**

For Research Use Only , Not for Diagnostic Use