


**Product Name: Recombinant Human HS1BP3 Protein**  
**Catalog #: PEH7740**



## Summary

<b>Name</b>	H1BP3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q53T59
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MQSPAVLVTS RRLQNAHTGL DLTVPQHQEY RGKMMSGHVE YQILVVTRLA AFKSAKHRPE DVVQFLVSKK YSEIEEFYQK LSSRYAAASL PPLPRKVLVY GESDIRERRA VFNEILRCVS KDAELAGSPE LLEFLGTRSP GAAGLTSRDS SVLDGTDSQT GNDEEAFDFF EEQDQVAEEG PPVQSLKGED AEESELEEEA LDPLGIMRSK KPKKHPKVAV KAKPSPRLTI FDEEVDPEDEG LFGPGRKLSP QDPSEDVSSV DPLKLFDDPD LGGAIPLGDS LLLPAACESG GPTPSLSHRD ASKELFRVEE DLDQILNLGA EPKPKPQLKP KPPVAAKPMI PRKPAVPPKA GPAEAVAGQQ KPQEQIQAMD EMDILQYIQD HDTPAQAAAPS LF

## image



## Background

**Alternative Names** HCLS1-binding protein 3 (HS1-binding protein 3) (HSP1BP-3)

**Product Name: Recombinant Human HS1BP3 Protein**  
**Catalog #: PEH7740**

---



**Note**

For Research Use Only , Not for Diagnostic Use