


Product Name: Recombinant Human CMPK2 Protein
Catalog #: PEH7741



Summary

Name	CMPK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5EBM0
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTRQALLELL GACQEAPRPH LGEFEADPRG QLWQRLWEVQ DGRRLQVGCA QVVPVPEPPL HPVVPDLPSS VVFPDREAAR AVLEECTSFI PEARAVLDLV DQCPKQIQKG KFQVVAIEGL DATGKTTVTQ SVADSLKAVL LKSPSCIGQ WRKIFDDEPT IIRRAFYSLG NYIVASEIAK ESAKSPVIDD RYWHSTATYA IATEVSGGLQ HLPPAHPVY QWPEDLLKPD LILLTVSPE ERLQRLQGRG MEKTREEAEL EANSVFRQKV

image



Background

Alternative Names UMP-CMP kinase 2, mitochondrial (EC 2.7.4.14) (Nucleoside-diphosphate kinase) (EC 2.7.4.6)

Product Name: Recombinant Human CMPK2 Protein
Catalog #: PEH7741



Note

For Research Use Only , Not for Diagnostic Use