


Product Name: Recombinant Human WDR44 Protein
Catalog #: PEH7742



Summary

Name	WDR44_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JSH3
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NGKAPDGQTV AGEVMGPQRP RSNSGRELTD EEILASVMIK NLDTGEEIPL SLAEEKLPTG INPLTLHIMR RTKEYVSNDA AQSDDDEKLQ SQPTDTDGGR LKQKTTQLKK FLGKSVKRAK HLAEEYGERA INKVKSVRDE VFHTDQDDPS SSDDEGMPYT RPVKFKAHAG FKGPDYFDQI KVVQDLSGEH MGAVWTMKFS HCGRLLASAG

image



Background

Alternative Names WD repeat-containing protein 44 (Rabphilin-11)

Note

Product Name: Recombinant Human WDR44 Protein
Catalog #: PEH7742



For Research Use Only , Not for Diagnostic Use