


Product Name: Recombinant Human CEP170 Protein
Catalog #: PEH7744



Summary

Name	CE170_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5SW79
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGGTRHRLPR EMIFVGRDDC ELMLQSRSDV KQHAVINYDA STDEHLVKDL GSLNGTFVND VRIPEQTYIT LKLEDKLRFG YDTNLFTVVQ GEMRVPEEAL KHEKFTIQLQ

image



Background

Alternative Names	Centrosomal protein of 170 kDa (Cep170) (KARP-1-binding protein) (KARP1-binding protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use