


Product Name: Recombinant Human FBNP1L Protein
Catalog #: PEH7745



Summary

Name	FBP1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T0N5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALEDFSHLPP EQRRKKLQQR IDELNRELQK ESDQKDALNK MKDVYEKNPQ MGDPGSLQPK LAETMNNIDR LRMEIHKNEA

image



Background

Alternative Names	Formin-binding protein 1-like (Transducer of Cdc42-dependent actin assembly protein 1) (Toca-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use