


Product Name: Recombinant Human SH3PXD2A Protein
Catalog #: PEH7750



Summary

Name	SPD2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TCZ1
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPKVRKKDA DANAGLTSCP RAKPSVRPKP FLNRAESQSQ EKMDISTLRR QLRPTGQLRG GLKGSKSEDS ELPPQTASEA

image



Background

Alternative Names	SH3 and PX domain-containing protein 2A (Adapter protein TKS5) (Five SH3 domain-containing protein) (SH3 multiple domains protein 1) (Tyrosine kinase substrate with five SH3 domains)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use