


**Product Name: Recombinant Human DDI2 Protein**  
**Catalog #: PEH7751**



## Summary

<b>Name</b>	DDI2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q5TDH0
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	45 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MLLTVYCVRR DLSEVTFSLQ VDADFELHNF RALCELESGI PAAESQIVYA ERPLTDNHRS LASYGLKDGD VVILRQKENA DPRPPVQFPN LPRIDFSSIA VPGTSSPRQR QPPGTQQSHS SPGEITSSPQ GLDNPALLRD MLLANPHEL LLKERNPPLA EALLSGDLEK FSRVLVEQQQ DRARREQERI RLFSADPFDL EAQAKIEEDI RQQNIEENMT IAMEEAPESF GQVVMLYINC KVNGHPVKAF VDSGAQMTIM SQACAERCNI MRLVDRRWAG IAKGVGTQKI IGRVHLAQQV IEGDFLPCSF SILEEQPMDM LLGLDMLKRH QCSIDLKKNV LVIGTTGSQT TFLPEGELPE CARLAYGAGR EDVRPEEIID QELAEALQKS AEDAERQKP

## image



## Background

**Alternative Names** Protein DDI1 homolog 2 (EC 3.4.23.-)

**Product Name: Recombinant Human DDI2 Protein**  
**Catalog #: PEH7751**

---



**Note**

For Research Use Only , Not for Diagnostic Use