


Product Name: Recombinant Human RAP1GAP2 Protein
Catalog #: PEH7754



Summary

Name	RPGP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q684P5
Host	E.coli
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YDEHEVNNTF KFGVIYQKAR QTLEEEELFGN NEESPAFKEF LDLLGDTITL QDFKGFRRGL DVTHGQTGVE SVYTTFRDRE IMFHVSTKLP FTDGDAQQLQ RKRHIGNDIV AIIFQEENTP FVPDMIASNF LHAYIVVQVE TPGTETPSYK VSVTAREDVP TFGPPLPSP VFQKGPEFRE FLLTKLTNAE NACCKSDKFA KLEDRTRAAAL

image



Background

Alternative Names Rap1 GTPase-activating protein 2 (Rap1GAP2) (GTPase-activating Rap/Ran-GAP domain-like protein 4)

Note

Product Name: Recombinant Human RAP1GAP2 Protein
Catalog #: PEH7754



For Research Use Only , Not for Diagnostic Use