


Product Name: Recombinant Human RAP1GAP2 Protein
Catalog #: PEH7755



Summary

Name	RPGP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q684P5
Host	E.coli
Species	Human
Molecular Mass	68 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEDKFENGGH GGFLESEFKRA IRVRSHSMET MVGGQKKSHS GGIPGSLSGG ISHNSMEVTK TTFSPPVVAA TVKNQSRSPI KRRSGLFPRL HTGSEGQGDS RARCDSTSST PKTPDGGHSS QEIKSETSSN PSSPEICPNK EKPFMKLKEN GRAISRSSSS TSSVSSTAGE GEAMEEGDSG GSQPSTTSPF KQEVFVYSPS PSSEPSLGA AATPIIMSRS PTDKSRNSP

image



Background

Alternative Names	Rap1 GTPase-activating protein 2 (Rap1GAP2) (GTPase-activating Rap/Ran-GAP domain-like protein 4)
--------------------------	---

Note

Product Name: Recombinant Human RAP1GAP2 Protein
Catalog #: PEH7755



For Research Use Only , Not for Diagnostic Use