


Product Name: Recombinant Human RAP1GAP2 Protein
Catalog #: PEH7756



Summary

Name	RPGP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q684P5
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YDEHEVNNTF KFGVIYQKAR QTLEELFGN NEESPAFKEF LDLLGDTITL QDFKGFRRGL DVTHGQTGVE SVYTTFRDRE IMFHVSTKLP FTDGDAQQLQ RKRHIGNDIV AIIFQEENTP FVPDMIASNF LHAYIVVQVE TPGTETPSYK VSVTAREDVP TFGPPLPSPP VFQKGPEFRE FLLTKLTNAE NACCKSDKFA KLEDRTRAAL

image



Background

Alternative Names	Rap1 GTPase-activating protein 2 (Rap1GAP2) (GTPase-activating Rap/Ran-GAP domain-like protein 4)
--------------------------	---

Note



Product Name: Recombinant Human RAP1GAP2 Protein
Catalog #: PEH7756

For Research Use Only , Not for Diagnostic Use