


Product Name: Recombinant Human STEAP4 Protein
Catalog #: PEH7757



Summary

Name	STEA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q687X5
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVCIFGTGDF GRSGLKMLQ CGYSVVFGR NPQKTTLLPS GAEVLSYSEA AKKSGIIIA IHREHYFLT ELTEVLNGKI LVDISNNLKI NQYPESNAEY LAHLVPGAHV VKAFNTISAW ALQSGALDAS RQVFVCGNDS KAKQRVMDIV RNLGLTPMDQ GSLMAAKEIE

image



Background

Alternative Names	Metalloreductase STEAP4 (EC 1.16.1.9) (Six-transmembrane epithelial antigen of prostate 4) (SixTransMembrane protein of prostate 2) (Tumor necrosis factor, alpha-induced protein 9)
--------------------------	--

Note

Product Name: Recombinant Human STEAP4 Protein
Catalog #: PEH7757



For Research Use Only , Not for Diagnostic Use