


Product Name: Recombinant Human ARHGAP17 Protein
Catalog #: PEH7758



Summary

Name	RHG17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68EM7
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LANQTVGRAE KTEVLSEDLL QIERRLDTVR SICHSHSKRL VACFQGQHG DAERRHKKLP LTALAQNMQE ASTQLEDSL L GKMLETCGDA ENQLALELSQ HEVFVEKEIV DPLYGIAEVE IPNIQQRKQ LARLVLDWDS VRARWNQAHK SSGTNFQGLP SKIDTLKEEM DEAGNKVEQC KDQLAADMYN FMAKEGEYGK FFVTLLEAQA DYHRKALAVL EKTLP EMRAH

image



Background

Alternative Names Rho GTPase-activating protein 17 (Rho-type GTPase-activating protein 17)
(RhoGAP interacting with CIP4 homologs protein 1) (RICH-1)

Note

Product Name: Recombinant Human ARHGAP17 Protein
Catalog #: PEH7758



For Research Use Only , Not for Diagnostic Use