


Product Name: Recombinant Human TWF2 Protein
Catalog #: PEH7760



Summary

Name	TWF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IBS0
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAHQTGIHAT EELKEFFAKA RAGSVRLIKV VIEDEQLVLG ASQEPVGRWD QDYDRAVLPL LDAQQPCYLL YRLDSQNAQG FEWLFLAWSP DNSPVRLKML YAATRATVKK EFGGGHIKDE LFGTVKDDLS FAGYQKHLSS CAAPAPL TSA ERELQQIRIN EVKTEISVES KHQTLQGLAF PLQPEAQRAL QQLKQKMOVNY IQMKDLERE TIELVHTEPT DVAQLPSRVP RDAARYHFFL YKHTHEGDPL ESVVFYISMP GYKCSIKERM LYSSCKSRL L DSVEQDFHLE IAKKIEIGDG AELTAEFLYD EVHPKQHAFK QAFAPKPGPG GKRGHKRLIR GPGENGDDS

image



Background

Alternative Names Twinfilin-2 (A6-related protein) (hA6RP) (Protein tyrosine kinase 9-like) (Twinfilin-1-like protein)

Product Name: Recombinant Human TWF2 Protein
Catalog #: PEH7760



Note

For Research Use Only , Not for Diagnostic Use